

# Cytokines

## Mouse Recombinant TRAIL



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

### TNF-related apoptosis-inducing ligand

Catalog # 78122  
78122.1

10 µg  
50 µg

## Product Description

TRAIL (TNF-related apoptosis-inducing ligand) belongs to the tumor necrosis factor (TNF) superfamily and is associated with initiating apoptosis. TRAIL has four major receptors: two death receptors DR4 and DR5, and two decoy receptors DcR1 and DcR2. TRAIL binds to the death receptors, which recruits the Fas-associated death domain (FADD) and activates caspases 8 and 10, which eventually leads to apoptosis (Pitti et al.; Wiley et al.; Zauli & Secchiero). It has been shown that mice lacking the expression of TRAIL have defects in thymocyte apoptosis and negative selection, and these mice had increased susceptibility to autoimmune diseases (Lamhamedi-Cherradi et al.).

## Product Information

**Alternative Names:** TL2, TNF-related apoptosis-inducing ligand, TNFSF10

**Accession Number:** P50592

**Amino Acid Sequence:** GRGGRPQKVA AHITGITRRS NSALIPISKD GKTLGQKIES WESSRKGHSF LNHVLFNRNGE LVIEQEGLYY IYSQTYFRFQ EAEDASKMVS KDKVRTKQLV QYIYKYTSYP DPIVLMKSAR NSCWSRDAEY GLYSIYQGGL FELKKNDRI F VSVTNEHLMD LDQEASFFGA FLIN

**Predicted Molecular Mass:** 20 kDa

**Species:** Mouse

**Cross Reactivity:** Human

**Formulation:** Lyophilized from a filtered solution in phosphate-buffered saline.

**Source:** E. coli

## Specifications

**Activity:** The specific activity is  $\geq 1.0 \times 10^4$  units/mg ( $EC_{50} \leq 100$  ng/mL) as determined by the ability to inhibit RPMI-8226 cell proliferation.

**Purity:**  $\geq 98\%$

**Endotoxin Level:** Measured by kinetic Limulus amebocyte lysate (LAL) analysis and is  $\leq 0.2$  EU/µg protein.


## Preparation and Storage

**Storage:** Store at  $-80^\circ\text{C}$ .

**Stability:** Stable as supplied for 12 months from date of receipt.

**Preparation:** Centrifuge vial before opening. Resuspend the product in sterile water or phosphate-buffered saline containing 0.1% bovine serum albumin (BSA) to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex. Store at  $2 - 8^\circ\text{C}$  for up to 1 week or at  $-20^\circ\text{C}$  to  $-80^\circ\text{C}$  for up to 3 months. Avoid repeated freeze-thaw cycles.  
NOTE: If reconstituted product will be used immediately BSA is not required.

## Data



(A) The biological activity of Mouse Recombinant TRAIL was tested by its ability to inhibit the proliferation of RPMI-8226 cells. Cell proliferation was measured using a fluorometric assay method. The EC<sub>50</sub> is defined as the effective concentration of the growth factor at which cell proliferation is at 50% of maximum. The EC<sub>50</sub> in the above example is less than 100 ng/mL.

(B) 2 µg of Mouse Recombinant TRAIL was resolved with SDS-PAGE under reducing (+) conditions and visualized by Coomassie Blue staining.

## Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, please visit our website at [www.stemcell.com/cytokines](http://www.stemcell.com/cytokines) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Lamhamedi-Cherradi S-E et al. (2003) Defective thymocyte apoptosis and accelerated autoimmune diseases in TRAIL<sup>-/-</sup> mice. *Nat Immunol* 4(3): 255–60.
- Pitti RM et al. (1996) Induction of apoptosis by Apo-2 ligand, a new member of the tumor necrosis factor cytokine family. *J Biol Chem* 271(22): 12687–90.
- Wiley SR et al. (1995) Identification and characterization of a new member of the TNF family that induces apoptosis. *Immunity* 3(6): 673–82.
- Zauli G & Secchiero P. (2006) The role of the TRAIL/TRAIL receptors system in hematopoiesis and endothelial cell biology. *Cytokine Growth Factor Rev* 17(4): 245–57.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.