

# Corning® 75 cm<sup>2</sup> Cell Culture, Straight Neck, Triangular Flask

**Polystyrene modified triangular, straight neck, cell culture flask with vent cap**

Catalog # 38014

100 Flasks



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

The Corning® 75 cm<sup>2</sup> Cell Culture, Straight Neck, Triangular Flask with vent cap is recommended for use with NeuroCult™ NS-A Proliferation Kit (Human) (Catalog #05751) and is also suitable for use with other cell culture media. This flask offers good pipette and cell scraper access to the corners, and the wide base provides added stability.

100 Flasks; 5 Flasks/bag (Corning Catalog #3276)

### APPLICATIONS

- Cell culture

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Flask: virgin polystyrene, hydrophilic and negatively charged
- Cap: high-density polyethylene with a 0.2 µm microporous polytetrafluoroethylene (PTFE) membrane
- Non-pyrogenic

### DIMENSIONS

- 155 x 84 x 40 mm (length x width x height)
- Good access to the corners
- Wide base provides added stability

### CAPACITY

- 90 mL
- Recommended working volume of 15 - 22.5 mL

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For information on chemical resistance or thermoplastic property, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and NeuroCult are trademarks of STEMCELL Technologies Canada Inc. Corning is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.