

# Costar® Microcentrifuge Tubes, 1.7 mL

## Polypropylene snap cap microcentrifuge tubes

Catalog # 38038

5000 Tubes



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Costar® Microcentrifuge Tubes are ideal for a wide variety of sample preparation, sample storage, and sample centrifugation tasks. The flat-top snap caps ensure a secure seal, are easy to open and close, and withstand repeated use. Tubes feature graduated markings and a frosted side panel and cap for easy labeling. Low-binding polymer reduces protein and nucleic acid binding to tube. Costar® Microcentrifuge Tubes are shipped non-sterile and are autoclavable.

0.65 mL - 500 Tubes; 500 Tubes/pack (Corning Catalog #3206 and equivalent to Corning Catalog #MCT-060-L-C)

1.7 mL - 5000 Tubes; 500 Tubes/pack (Corning Catalog #3621 and equivalent to Corning Catalog #MCT-175-L-C)

### APPLICATIONS

- Liquid handling
- PCR
- Molecular biology
- Cell culture
- FACS

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Tube and cap: clear and transparent polypropylene
- Easy-open snap cap
- RNase/DNase-free and non-pyrogenic

### DIMENSIONS

- 12.7 x 9.2 x 39.0 mm (outer diameter x inner diameter x height)

### NOMINAL VOLUME

- 1.7 mL
- Graduated markings

### CENTRIFUGAL FORCE

- Up to 17000 x g

### STERILITY

- Non-sterile; autoclavable

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Costar is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.